



# Mixing Active Pharmaceutical Ingredients (APIs):

Process, Risks, Cleaning,  
and How PerMix Solves It


# Why API Mixing Is Harder Than It "Should" Be

APIs behave like the universe's way of reminding humans that small things can be terrifying.

Most pharmaceutical blends combine:

- ‡ Micronized APIs (often cohesive, electrostatic, poorly flowing)
- ‡ Excipients (often larger, denser, free-flowing)
- ‡ Low inclusion rates (sometimes <1%, sometimes ppm-level with potent compounds)

## That combination triggers four classic failure mechanisms:



### Segregation (the silent assassin)

Even if you achieve a great blend inside the mixer, segregation can happen:

- ‡ During transfer (vacuum conveying, gravity drops)
- ‡ During discharge (funnel flow vs mass flow)
- ‡ During bin filling (particle percolation and air entrainment)
- ‡ During vibration (equipment, forklifts, packaging machines)



### Electrostatics (powder cling + dose drift)

Static causes:

- ‡ API coating on walls, shafts, seals
- ‡ inconsistent discharge weights
- ‡ "hidden hold-up" that contaminates the next batch or steals yield



### Agglomeration (fake uniformity)

Fine APIs can "look" blended while remaining as micro-clusters. Those clusters become content-uniformity failures when sampled.



### Over-shear / under-shear (Goldilocks energy problem)

Too much energy: attrition, heat, PSD shift, possible polymorph risk for sensitive materials.

Too little energy: hot spots, agglomerates survive, CV doesn't converge.

# The API Mixing Workflow (Done Properly)

## Step

1

### Raw Material Conditioning

Before mixing even begins, many plants:

- ▄ de-lump APIs (sieve or gentle deagglomeration)
- ▄ control humidity (static control + flow consistency)
- ▄ manage temperature (some APIs are moisture/heat sensitive)

PerMix approach: choosing mixer type and options based on whether the powder is cohesive, segregative, fragile, or abrasive—because "powder" is not one material. It's a behavior.

## Step

2

### Charging Strategy

**How You Load Matters as Much as How You Mix**

Bad charging causes immediate stratification:

- ▄ dumping API on top → it floats and smears
- ▄ adding API early into high-flow excipients → it gets buried and never disperses
- ▄ vacuum loading without controls → API fines deposit in lines and filters

Best-practice charging methods include:

- ▄ geometric dilution (API first blended into a "pre-blend" excipient fraction)
- ▄ split charging (API introduced in multiple additions)
- ▄ contained addition (isolator, glovebox interface, split butterfly valve)
- ▄ vacuum loading with dust-safe filtration and controlled inlet design

PerMix can support these strategies through inlet placement, loading accessories, and designs that reduce dead zones and deposition.

## Step

3

### Mixing Phase

#### Particle Engagement, Not Just Bulk Circulation

This is the core truth: API uniformity requires particle-level interaction, not just moving a pile.

#### Why legacy mixers struggle

- ▄ **Ribbon mixers** are powerful but can create preferential flow paths and dead zones—especially with cohesive APIs.
- ▄ **Tumble blenders** are gentle but often slow and sensitive to fill level and PSD mismatch; liquid addition can be tricky.

#### How PerMix answers with the right mixing physics

##### Fluidized Zone Mixing

When the formulation tends to segregate or includes low-dose actives.

Instead of relying on gravity circulation alone, fluidized-zone designs keep the powder bed continuously reoriented. This helps:

- ▄ keep light and heavy fractions engaged
- ▄ reduce density-driven stratification
- ▄ break weak agglomerates gently
- ▄ converge CV faster and more reliably

##### Paddle / Plow Mixing with Intensification

When you need controlled shear + fast convergence.

For blends that need deagglomeration without overworking the entire batch, PerMix integrates:

- ▄ choppers / intensifiers positioned to target clusters
- ▄ adjustable speeds to dial energy input
- ▄ options for liquid dispersion where required

## Step

4

### Liquid Addition (If Applicable)

#### The "Wet Spot" Trap

Some pharma processes require binders, granulation liquids, or functional coating solutions. The #1 mistake is adding liquid as a stream.

Proper liquid addition requires:

- ▄ atomization (fine droplets)
- ▄ correct nozzle placement (into active mixing zone, not a dead pocket)
- ▄ rate control (avoid overwetting and secondary agglomeration)
- ▄ synchronization (liquid addition + mixing energy profile)

PerMix designs can incorporate spray manifolds and control-friendly systems so "liquid addition" becomes engineered—not improvised.

## Step

5

### Discharge and Transfer

#### Where Great Blends Go to Die

Even perfect in-mixer uniformity can be destroyed by discharge.

Key discharge risks:

- ▄ ratholing and funnel flow
- ▄ segregation in downspouts
- ▄ sifting segregation during vibration
- ▄ product hold-up that contaminates the next batch

PerMix design focus:

- ▄ discharge geometry that supports consistent flow
- ▄ options that reduce hold-up zones
- ▄ cleaning-aware outlet designs (because discharge is also a cleaning hotspot)

# Validation, Sampling, and Why “It Looked Fine” Isn’t a Spec

## API mixing performance is measured with:

- ▄ content uniformity testing
- ▄ blend uniformity sampling plans
- ▄ CV targets (especially for low-dose actives)
- ▄ sometimes PAT tools (process analytical technology) depending on plant strategy

## But equipment matters because:

- ▄ dead zones create hidden failure modes
- ▄ inconsistent hold-up breaks repeatability
- ▄ poor cleanability increases cross-contamination risk and downtime

Which brings us to the part most plants feel in their bones: cleaning.

## Cleaning Pharmaceutical Mixers: Wet CIP, Dry CIP, and Hybrid Strategies

Cleaning is not "maintenance." In GMP pharma, cleaning is a validated manufacturing step with compliance implications.

PerMix supports multiple cleaning philosophies depending on:

- ▄ product toxicity/potency
- ▄ solubility characteristics
- ▄ allergen/cross-contamination requirements
- ▄ available utilities (water, steam, air, vacuum)
- ▄ turnaround time expectations



# Wet CIP (Clean-in-Place):

## When You Need Washdown Certainty

Wet CIP uses water (often with detergent) delivered through spray devices to remove residues.

### Where wet CIP shines

- ▣ sticky residues
- ▣ water-soluble products
- ▣ products that require sanitization
- ▣ facilities standardized on wash validation methods

### What makes wet CIP succeed

Wet CIP is only as good as coverage + drainability. A "CIP system" that can't hit shadow areas is just a confidence machine.

PerMix CIP engineering typically focuses on:

- ▣ spray device placement and coverage of the full internal envelope
- ▣ cleaning access to shafts, seals, outlet, and intensifier housings
- ▣ drain design to prevent pooling
- ▣ surface finishes that reduce adhesion and speed rinse

### Real-world constraints

Wet CIP consumes:

- ▣ water
- ▣ effluent handling capacity
- ▣ time for drying (which becomes a microbial/validation concern in some environments)

So wet CIP is powerful—but not always the fastest path back to production.



## Dry CIP:

## When Downtime and Utilities Matter

Dry CIP uses directed air (and/or vacuum-assisted) systems to dislodge and evacuate powder residues without water.

### Where dry CIP shines

- ▣ dry powders that are non-sticky
- ▣ rapid product changeovers
- ▣ plants trying to reduce wastewater and downtime
- ▣ products that clump or degrade with moisture
- ▣ operations where "wet cleaning" creates drying bottlenecks

PerMix dry CIP advances (the practical version):

- ▣ strategically placed air nozzles to hit typical build-up zones
- ▣ air + vacuum logic to remove loosened fines rather than just re-suspending them
- ▣ designs that avoid "dust traps" and ledges

Dry CIP isn't just a blower. The trick is to move residues out of the machine, not into the air.

# Hybrid Cleaning (Dry + Wet):

## The Fastest Path to Control

This is often the best of both worlds.

A hybrid strategy typically works like:

1. **Dry CIP** first to remove bulk powder hold-up quickly (reduces what wet CIP must dissolve)
2. **Wet CIP** second for final residue removal, sanitization, and validated endpoint
3. Optional drying/air purge to return the mixer to service faster

### Benefits:

- ✦ less water usage than full wet CIP alone
- ✦ faster turnaround than wet-only cleaning
- ✦ reduced effluent load
- ✦ better control over residues that are partly soluble and partly adherent

In practice, hybrid cleaning can be a major productivity lever, especially in multi-SKU operations.

## Putting It Together:

## How PerMix Answers the Full API Mixing + Cleaning Problem

### Outcome targets:

- ✦ reliable low-dose homogeneity
- ✦ controlled energy input (no overworking)
- ✦ repeatable discharge behavior
- ✦ minimized hold-up and dusting
- ✦ cleanability engineered into the geometry
- ✦ CIP strategies aligned with your plant constraints (wet, dry, hybrid)
- ✦ reduced downtime between SKUs

Because the real KPI isn't "does it mix."

The KPI is validated batches per week, with minimal deviation risk, and fast, documented changeovers.

## PerMix Easy-Clean Mixer Options Designed for Reality


Even the best CIP system can't compensate for a mixer full of crevices.

Easy-Clean is about removing the places residue likes to hide:

- ✦ simplified internal geometry
- ✦ reduced ledges and shadow zones
- ✦ clean-friendly shaft and seal arrangements
- ✦ faster access to internal tools (intensifiers/choppers) depending on configuration
- ✦ surfaces and weld practices aimed at cleanability and repeatability

In a pharma setting, Easy-Clean design translates to:

- ✦ less manual intervention
- ✦ shorter cleaning cycles
- ✦ more reliable verification (visual + swab)
- ✦ reduced cross-contamination risk
- ✦ higher uptime and more predictable scheduling



## Global Factory & Office

**PerMix North America (Factory)**  
Chicago, Illinois, 60601, USA  
Contact: John Paul  
Tel: +1 630-649-1357  
Email: [john.paul@permixtec.com](mailto:john.paul@permixtec.com)  
[www.permixmixers.com](http://www.permixmixers.com)

**PerMix Europe (Factory)**  
9900 Eeklo - Belgium  
Contact: Dimitri Baeten  
Tel: +32 9 430 50 17  
Email: [dimitri.eu@permixtec.com](mailto:dimitri.eu@permixtec.com)  
[www.permixmixers-eu.com](http://www.permixmixers-eu.com)

**PerMix China (Factory)**  
Shanghai, 201821, China  
Contact: Aaron Huang  
Tel: +86 133 911 68218  
Email: [aaron.h@permixtec.com](mailto:aaron.h@permixtec.com)  
[www.permix.cn](http://www.permix.cn)

**PerMix India**  
Navi Mumbai, 400 708, India  
Contact: Mr. Anant Wagh  
Tel: +91 99304 69228  
Email: [india@permixtec.com](mailto:india@permixtec.com)  
[www.permixmixers.in](http://www.permixmixers.in)

**PerMix UK**  
Scotland, Great Britain  
Contact: Laura Griffin  
Tel: +44 7554 139667  
Email: [laura.griffin@permixtec.com](mailto:laura.griffin@permixtec.com)  
[www.permixmixers.co.uk](http://www.permixmixers.co.uk)

**PerMix Ukraine**  
Kyiv, 02095, Ukraine  
Contact: Viacheslav Pishyi  
Tel: +38 067 45 44 694  
Email: [ua@permixtec.com](mailto:ua@permixtec.com)  
[www.permixtec.com](http://www.permixtec.com)

## Headquarter

**PerMix Israel**  
Adi 17940, Israel  
Contact: Mr. Arie Srugo  
Tel: +972-54-908-0144  
Email: [srugo@permixtec.com](mailto:srugo@permixtec.com)  
[www.permixtec.com](http://www.permixtec.com)